


## Ferromagnetic Resonance (FMR) & Inverse Spin Hall Effect (ISHE) Insert

The Ferromagnetic Resonance (FMR) probe is designed to accept a Coplanar Waveguide that allows the measurement of in-plane and out-of-plane ferromagnetic resonance spectroscopy using the same insert. The unit is supplied with Cryogenic semi-rigid coaxial cable assemblies, associated RF connectors and an easy and repetitive sample mounting fixture.



Frequency Sweep at Room Temperature, showing result at 5 GHz as an example

### Key Features:

- » **Dual FMR Probes:**  
Two insert FMR probes with co-planar waveguides (CPW) for in-plane and out-of-plane FMR and in-plane ISHE measurements
- » **ISHE-Ready:**  
DC lines for ISHE measurements, with precision alignment stage
- » **Modulation Coils:**  
Surround the CPW and sample platform, enabling standard FMR with magnetic field modulation and Schottky diode detection
- » **Broadband Capability:**  
Perform broadband measurements with a vector network analyser (VNA)
- » **Dedicated Software:**  
To perform frequency dependent FMR measurement series (i.e. resonance field vs. frequency and field line width vs. frequency).



## Unlock precision and versatility in spintronics research

### Precision Sample Alignment for ISHE

- » Direct electrical contacting of the sample for inverse spin Hall effect (ISHE) measurements through a specialised alignment stage;
- » Enables accurate sample alignment relative to the external magnetic field;
- » An optical microscope combined with a linear translation stage ensures precise positioning of the device under test (DUT);
- » Achieve angular precision better than 1° between the DUT and CPW.

### Performance:

Bandwidth	Up to 40 GHz
Microwave Signal Generator	To meet research requirements
Lock-in Amplifier	MFLI or similar (optional upgrade to VNA)
Magnetic Field	Up to ±16 Tesla
Temperature Range	4 - 300 K
Magnetic Parameters	<ul style="list-style-type: none"> <li>• Field linewidth</li> <li>• Damping constant</li> <li>• Sample saturation magnetization</li> <li>• Anisotropy constant</li> </ul>
Sample Size	Up to 10 mm <sup>2</sup>
DC Lines for ISHE Measurements	2+4

