

# CRYOGEN FREE ESR / EPR MEASUREMENT SYSTEM



- **Magnetic fields to 16 Tesla**
- **High homogeneity 5-ppm magnet**
- **Optional sweep coil giving +/- 0.006T (optional)**
- **Fitted with leads and a switch for persistent mode operation.**
- **Integrated Variable Temperature Insert (optional)**
- **Low operating costs - no liquid helium required**
- **User-friendly and convenient operation**
- **Fully automated**

CRYOGENIC Ltd.  
30 Acton Park Industrial Estate  
The Vale, London W3 7QE, UK  
Tel: +44 (20) 8743 6049  
Fax: +44 (20) 8749 5315  
E-mail: [info@cryogenic.co.uk](mailto:info@cryogenic.co.uk)

[www.cryogenic.co.uk](http://www.cryogenic.co.uk)

**CRYOGENIC**  
CRYOGENIC LIMITED

[www.cryogenic.co.uk](http://www.cryogenic.co.uk)